

Analog Dual UST

an unique System for Physiotherapy Applications



FEATURES

- Feather touch keypad
- Ruggedized probe holder
- SMPS powered
- Single probe for 1 and 3 MHz output
- Configurable treatment cycle and Intensity level
- Stand alone START/STOP switch for output control
- User friendly selector knobs for therapy mode selection
- User friendly 7 segment display for treatment cycle and output power indication
- Equipped with "PIC" microcontroller

SPECIFICATIONS

- Output frequency : 1 and 3MHz
- Output Mode : Continuous & Pulse
- Pulse Ratio : 1:1,1:2, 1:4, 1:7,1:8 & 1:16
- Power Level : 0.1 to 3.1 watts/cm²
- Timer : 1 – 15 minutes adjustable for treatment time control
- Input Power : SMPS @ 100V - 275V AC, 50 Hz Operated
- Dimension : 185 x 205 x 90 mm
- Weight : 1.6 kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics