

Analog Traction

an unique System for Physiotherapy Applications



FEATURES

- Configurable treatment cycle
- Emergency Stop control for patient safety
- Stand alone START/STOP switch for control
- User friendly selector knob for therapy selection
- LED for treatment cycle and output indication
- Equipped with "PIC" microcontroller

SPECIFICATIONS

- Therapy Modes : Lumbar & Cervical
- Weight : Lumbar 20 to 40Kg at 2Kg scale
Cervical 4 to 14Kg at 1 kg scale
- Hold time : 10,20,40,60,80 Sec
- Rest time : 1,5,10,15,20 Sec
- LED Indicator : Available for "Hold and Rest" time
- Timer : 1 – 30 minutes adjustable
- Power Supply : LPS @ 230VAC, 50Hz
- Dimension : 180 x 260 x 260 mm
- Weight : 10 Kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics