

DIGITAL LASER THERAPY

an unique System for Physiotherapy Applications



FEATURES

- Feather-touch keypad
- Ruggedized probe holder
- 75 therapeutic-library programs
- Graphical LCD display (128 x 64)
- SMPS powered
- Equipped with "PIC" microcontroller
- Single therapy probe for IR and RED
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Energy level are configurable

SPECIFICATIONS

- IR wavelength : 780 nm
- RED wavelength : 660 nm
- Energy Level : 0.1 to 25 mJ
- Timer : 1 to 8 min adjustable
- Power Supply : SMPS @ 100V-275VAC, 50Hz
- Dimension : 250 x 250 x 110mm
- Weight : 2.17 kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics