

Digital Traction

an unique System for Physiotherapy Applications



FEATURES

- ▶ Feather-touch keypad
- ▶ Graphical LCD display (128 x 64)
- ▶ Treatment cycle is configurable
- ▶ Equipped with "PIC" microcontroller
- ▶ 20 configurable Favorites memory location
- ▶ Emergency Stop switch for Patient safety
- ▶ Stand alone Start/ Stop switch for output control

SPECIFICATIONS

- ▶ Therapy Modes : Static & Intermittent
- ▶ Weight : Cervical: 4 to 14 Kg at 1 Kg scale
Lumbar: 20 to 40 Kg at 2 Kg scale
- ▶ Hold time : 5 to 99 Sec @ 5 Sec scale
- ▶ Rest time : 1 to 99 Sec @ 5 Sec scale
- ▶ Timer : 1 to 30 min adjustable
- ▶ Power Supply : LPS @ 230VAC, 50Hz
- ▶ Dimension : 280 x 260 x 260 mm
- ▶ Weight : 10 Kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics