

# STIMULATOR - T

an unique System for Physiotherapy Applications



## FEATURES

- Configurable Intensity level
- Independent LED for output indication
- User friendly selector knob for therapy selection
- Equipped with "PIC" microcontroller
- Buffer indication for IG-ON Time

## SPECIFICATIONS

- Therapy Modes : Faradic, Surge Faradic, Galvanic, Interrupted Galvanic
- LED Indicator : for Faradic & Surge Faradic, Galvanic & Interrupted Galvanic
- Faradic : 0.7mS ON time @20mS period
- Surge Faradic : Faradic @ Hold & Rest time:1 to 6 sec adjustable
- Interrupted Galvanic : ON time : 0.01, 0.03, 0.1, 0.3, 1,3, 10, 30, 100, 300&1000 ms  
Rest time : 1 to 6 sec adjustable
- Output : 0 – 99 mA adjustable
- Power Supply : LPS @ 230VAC, 50Hz
- Dimension : 185 x 205 x 90 mm
- Weight : 2.20 Kg



No.6/1, Venkatachalam Street, Pari Nagar,  
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com  
www.alphameditronics.in



ALPHA

## ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147  
+91-73586 20247



alphameditronics