

# Analog 4 Channel TENS

an unique System for Physiotherapy Applications



## FEATURES

- Four output channels with isolated output
- Configurable frequency and intensity level for each channel
- Configurable HOLD and REST time
- User friendly selector switch for therapy selection
- Stand alone START/STOP switch for control
- Feather touch keypad for treatment cycle
- LED for frequency and output indication

## SPECIFICATIONS

- Mode : Continuous & Burst
- Frequency : up to 150 Hz @ 400  $\mu$ s ON time
- Burst period : HOLD & REST: 1 to 6 Sec adjustable
- Output : 0 – 65mA adjustable
- Timer : 1 – 30 minutes adjustable
- Power Supply : SMPS @ 100V–275V AC, 50Hz
- Dimension : 280 x 170 x 80 mm
- Weight : 1.5 Kg



No.6/1, Venkatachalam Street, Pari Nagar,  
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com  
www.alphameditronics.in



ALPHA

**ALPHA MEDITRONICS**

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147  
+91-73586 20247



alphameditronics