

Analog mini IFT

an unique System for Physiotherapy Applications



SPECIFICATIONS

- Medium Frequency : 2, 3, 4, 5, 6 KHz and Russian Current
- Base : up to 150Hz
- Spectrum : up to 100Hz with Linear, Non-Linear and Trapezoidal
- Mode : 4P, 4PV, 2P
- Vector type : Linear, Non-Linear & Trapezoidal
- Output : Isolated 0 – 100V Output with display meter & LED Indicator
- Timer : 1 – 30 minutes adjustable
- Input Power : SMPS @ 100V – 275V AC, 50 Hz Operated
- Dimension : 250 x 235 x 75 mm
- Weight : 2.0 kg

FEATURES

- Feather touch keypad
- SMPS powered
- Isolated dual output channels
- Configurable treatment cycle and Intensity level
- Stand alone Start/ Stop switch for output control
- User friendly selector knobs for therapy mode selection
- User friendly 7 segment display for treatment cycle and output indication
- Equipped with "ARM" microcontroller



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics