

mini UST

an unique System for Physiotherapy Applications



FEATURES

- Feather touch keypad
- SMPS powered
- Isolated probe with elevated safety
- Ruggedized probe holder
- Configurable treatment cycle and Intensity level
- Stand alone Start/ Stop switch for output control
- User friendly selector knob for therapy mode selection
- User friendly 7 segment display for treatment cycle and output indication
- Equipped with "PIC" microcontroller

SPECIFICATIONS

- Output frequency : 1MHz
- Output Mode : Continuous & Pulse
- Pulse Ratio : 1:1,1:2, 1:4, 1:7,1:8 & 1:16
- Power Level : 0.1 to 3.1 watts/cm²
- Timer : 1 – 15 minutes adjustable
- Power Supply : SMPS @ 100V-275VAC, 50Hz
- Dimension : 165 x 165 x 75mm
- Weight : 1.4 kg

No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083

+91-44-4775 9102

+91-73586 20147 / 73586 20247

alphameditronics@gmail.com

www.alphameditronics.com
www.alphameditronics.in



ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics