

Digital 4 Channel TENS

an unique System for Physiotherapy Applications



FEATURES

- Feather-touch keypad
- Graphical LCD display (128 x 64)
- SMPS powered
- Equipped with "PIC" microcontroller
- 4 independent channels with isolated output
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Intensity level are configurable
- User friendly knob to set intensity

SPECIFICATIONS

- Mode : Continuous & Burst
- Frequency : up to 150 Hz @ 400 μ s ON time
- Burst period : HOLD & REST: 1 to 6 Sec adjustable
- Output : 0 – 65mA adjustable
- Timer : 1 – 30 minutes adjustable
- Power Supply : SMPS @ 100V–275V AC, 50Hz
- Dimension : 185 x 205 x 90 mm
- Weight : 1.75 Kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics