

DIGITAL IFT

an unique System for Physiotherapy Applications



FEATURES

- Feather-touch keypad
- 100 therapeutic-library programs
- Graphical LCD display (128 x 64)
- SMPS powered
- Isolated probes with elevated safety
- Stand alone Start/ Stop switch for control
- 50 configurable Favorites memory location
- User friendly knob to set intensity
- Equipped with "PIC" microcontroller

SPECIFICATIONS

- Frequency : 2, 3, 4, 5 and 6 KHz
- Base : 1 to 150Hz @ 1 Hz step
- Spectrum : 1 to 150Hz @ 1 Hz step with Linear, Non-Linear and Trapezoidal shape
- Mode : 4P, 4PV, 2P and Russian Current
- Vector type : 2PL, 4PL, 4PT & 4PNL
- Output : Isolated 0 – 100V
- Timer : 1 – 30 minutes adjustable
- Input Power : SMPS @ 100V – 275V AC, 50 Hz Operated
- Dimension : 185 x 205 x 90 mm
- Weight : 2 kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics