

DIGITAL MST-TENS

an unique System for Physiotherapy Applications



FEATURES

- Feather-touch keypad, SMPS powered
- Equipped with "PIC" microcontroller
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Output level are configurable
- 25 therapeutic-library programs
- 20 configurable Favorites memory locations
- Big 128 x 64 Graphical LCD display
- User friendly knob to set intensity

SPECIFICATIONS

- Therapy Modes : Faradic, Surge Faradic, Galvanic, Interrupted Galvanic, U- Reiz, HVPC, Micro Current & TENS
- TENS
 - » Mode : Continuous & Burst
 - » Frequency : 1 to 150 Hz @ 25 to 400 μ s ON time adjustable
 - » Burst period : Hold & Rest @ 1-6 sec adjustable
 - » Output : 0 – 100 mA adjustable
- STIMULATOR
 - » Faradic : 0.7mS ON time @20mS period
 - » Surge Faradic : Faradic @ Hold & Rest :1 to 6 sec adjustable
 - » Interrupted Galvanic : ON time : 0.01,.03,0.1,0.3,1,3,10,30,100,300&1000 ms
Rest time :1-6 sec adjustable
 - » Output Intensity : 0 – 100mA adjustable
- Timer : 1 to 30 min adjustable
- Power Supply : SMPS @ 100V–275VAC, 50Hz
- Dimension : 185 x 205 x 90 mm
- Weight : 1.75 Kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics