

# DIGITAL TENS Plus

an unique System for Physiotherapy Applications



## FEATURES

- Isolated Dual output channels
- Configurable frequency and intensity level using 7 segment display and Feather-touch Keys
- Equipped with "PIC" microcontroller
- Stand alone Start/ Stop switch for output control
- Built-in 15 minutes timer

## SPECIFICATIONS

- Mode : TENS: Continuous, Burst and 5 sweep modes  
MST: Faradic, Surge Faradic, Galvanic and Interrupted Galvanic
- Frequency : 1 to 150 Hz with 5 Hz step adjustable @ 400  $\mu$ s ON time
- Sweep frequency : 1 to 50Hz, 50 to 100Hz, 100 to 150Hz, 1 to 100Hz and 1 to 150Hz with Continuous and Burst mode.
- Burst period : Hold @ 3 sec, Rest @ 3 Sec
- Output : 0 – 65 mA adjustable
- Power Supply : 12V Adaptor Operated
- Dimension : 95 x 115 x 35 mm
- Weight : 300 g



No.6/1, Venkatachalam Street, Pari Nagar,  
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com  
www.alphameditronics.in



ALPHA

**ALPHA MEDITRONICS**

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147  
+91-73586 20247



alphameditronics