


DIGITAL UST+

an unique System for Physiotherapy Applications



FEATURES

- Feather-touch keypad
- Ruggedized Probe holder
- 75 therapeutic-library programs
- 20 configurable Favorites memory locations
- Big 128 x 64 Graphical LCD display
- Isolated probe with elevated safety
- SMPS powered
- Equipped with "PIC" microcontroller
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Intensity level are configurable

SPECIFICATIONS

- Output frequency : 1MHz
- Output Mode : Continuous & Pulse
- Pulse Ratio : 16 modes from 1:1 to 1:16
- Power Level : 0.1 to 3.1 watts/cm²
- Timer : 1 – 30 minutes adjustable
- Power Supply : SMPS @ 100V–275VAC, 50Hz
- Dimension : 185 x 205 x 90 mm
- Weight : 1.9 kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics