

STIMULATOR Plus

an unique System for Physiotherapy Applications



FEATURES

- Feather touch keypad
- SMPS powered
- Stand alone Start/ Stop switch for output control
- Configurable treatment cycle and Intensity level
- User friendly selector knobs for therapy mode selection
- User friendly 7 segment display for treatment cycle and output indication
- Equipped with "PIC" microcontroller

SPECIFICATIONS

- Therapy Modes : F, SF, G, IG, U-REIZ, HVPC, Micro Current
- Faradic : 0.7mS ON time @20mS period
- Surge Faradic : Faradic @ Hold & Rest:1 to 6 sec adjustable
- Interrupted Galvanic : ON time : 0.1, 0.3, 1, 3, 10, 30, 100, 300 & 1000 ms
Rest time : 1 to 6 sec adjustable
- Output : 0 – 65 mA adjustable
- Timer : 1 – 30 minutes adjustable
- Power Supply : SMPS @ 100V–275VAC, 50Hz
- Dimension : 185 x 205 x 90 mm
- Weight : 1.7Kg



No.6/1, Venkatachalam Street, Pari Nagar,
Jafferkhanpet, Chennai- 600083



+91-44-4775 9102



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics