

# TENS Combo

an unique System for Physiotherapy Applications



## FEATURES

- Isolated Dual output channels
- Configurable frequency and intensity level
- User friendly selector switch for therapy selection
- LED for frequency and output indication
- TENS, Faradic and Interrupted Galvanic modes in single handy machine.

## SPECIFICATIONS

- Mode : Continuous & Burst
- Frequency : up to 150 Hz @ 400  $\mu$ s ON time
- Burst period : Hold @ 3 sec, Rest @ 3 Sec
- Output : 0 – 65 mA adjustable
- Power Supply : 9V battery / Adaptor Operated
- Dimension : 105 x 100 x 60 mm
- Weight : 400 g

No.6/1, Venkatachalam Street, Pari Nagar,  
Jafferkhanpet, Chennai- 600083

+91-44-4775 9102

+91-73586 20147 / 73586 20247

alphameditronics@gmail.com

www.alphameditronics.com  
www.alphameditronics.in



## ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147  
+91-73586 20247



alphameditronics