

# Wax Bath Therapy

an unique System for Physiotherapy Applications



## CAPACITY: 1 KG

- Weight: 3.75 Kg
- Dimension: 24 x 24 x 16.5 Cm



## CAPACITY: 3 KG

- Weight: 5.10 Kg
- Dimension: 26 x 26 x 19 Cm




## CAPACITY: 5 KG

- Weight: 6.70 Kg
- Dimension: 35.5 x 28 x 19 Cm

### Features and Specifications:


- MS body with wooden top
- Comes with ON-OFF switch
- Indicator LED for mains and Thermostat ON
- Fuse for protection
- Thermostat control up to 110 C
- Operating voltage: 230 V AC, 50 Hz

 No.6/1, Venkatachalam Street, Pari Nagar,  
Jafferkhanpet, Chennai- 600083

 +91-44-4775 9102

 +91-73586 20147 / 73586 20247

 alphameditronics@gmail.com

 www.alphameditronics.com  
www.alphameditronics.in



ALPHA

## ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147  
+91-73586 20247



alphameditronics