

Semi Digi IFT

an unique System for Physiotherapy Applications



FEATURES

- Isolated dual channel Output
- SMPS powered
- Equipped with "ARM" microcontroller
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Intensity level are configurable
- User friendly 7 segment display, Bar Graph LED and Feather-touch Keys as user interface

SPECIFICATIONS

- Medium Frequency : 2, 3, 4, 5, 6 KHz and Russian
- Base : up to 150Hz
- Spectrum : up to 100Hz with Linear, Non-Linear and Trapezoidal shape
- Mode : 4P, 4PV, 2P
- Vector type : Linear, Non-Linear & Trapezoidal
- Output : Isolated 0 – 100V
- Timer : 1 – 30 minutes adjustable
- Input Power : SMPS @ 100V – 275V AC, 50 Hz Operated
- Dimension & Weight : 290 x 290 x 110 mm, 2.6 Kg



New No.8, Old No. J71, Karikalan Street
Jafferkhanpet, Chennai - 600 083



+91-44-4858 2421



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics