

Semi Digi Laser

an unique System for Physiotherapy Applications



FEATURES

- Ruggedized probe holder
- SMPS powered
- Equipped with "PIC" microcontroller
- Single therapy probe for IR and RED
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Energy level are configurable
- User friendly 7 segment display, Bar Graph LED and Feather-touch Keys as user interface

SPECIFICATIONS

- IR wavelength : 780 nm
- RED wavelength : 660 nm
- Energy Level : 0.1 to 25 mJ
- Power Supply : SMPS @ 150V-275VAC, 50Hz
- Timer : 1 to 8 min adjustable
- Dimension : 290 x 290 x 160mm
- Weight : 2.58 kg



New No.8, Old No. J71, Karikalan Street
Jafferkhanpet, Chennai - 600 083



+91-44-4858 2421



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics