

Semi Digi Stim

an unique System for Physiotherapy Applications



FEATURES

- 9-therapeutic currents offered
- SMPS powered
- Equipped with "PIC" microcontroller
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Intensity level are configurable
- User friendly 7 segment display, Bar Graph LED and Feather-touch Keys as user interface

SPECIFICATIONS

- Therapy Modes : Faradic, Surge Faradic, Galvanic, Interrupted Galvanic, U- Reiz, HVPC, Micro Current, & TENS
- TENS
 - » Mode : Continuous & Burst
 - » Frequency : up to 150 Hz @ 400 μ s ON time
 - » Burst period : Hold & Rest @ 1-6 sec adjustable
 - » Output : 0 – 65mA adjustable and isolated
- STIMULATOR
 - » Faradic : @ Hold & Rest: 1 to 6 sec adjustable
 - » Surge Faradic : Faradic @ Hold & Rest: 1 to 6 sec adjustable
 - » Interrupted Galvanic : ON time : 0.01, .03, 0.1, 0.3, 1, 3, 10, 30, 100, 300 & 1000 ms
Rest time : 1-6 sec adjustable
 - » Output Intensity : 0 – 65mA adjustable
- Timer : 1 to 30 min adjustable
- Power Supply : SMPS @ 150V-275VAC, 50Hz
- Dimension : 290 x 260 x 160 mm
- Weight : 2.7 Kg



New No.8, Old No. J71, Karikalan Street
Jafferkhanpet, Chennai - 600 083



+91-44-4858 2421



+91-73586 20147 / 73586 20247



alphameditronics@gmail.com



www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics