

Semi Digi Ultrasound

an unique System for Physiotherapy Applications




FEATURES


- Ruggedized probe holder
- Isolated probe with elevated safety
- SMPS powered
- Equipped with "PIC" microcontroller
- Stand alone Start/ Stop switch for output control
- Treatment cycle and Intensity level are configurable
- User friendly 7 segment display, Bar Graph LED and Feather-touch Keys as user interface

SPECIFICATIONS


- Output frequency : 1MHz
- Output Mode : Continuous & Pulse
- Pulse Ratio : 1:1,1:2, 1:4, 1:7,1:8 & 1:16
- Power Level : 0.1 to 3.1 watts/cm²
- Power Supply : SMPS @ 150V-275VAC, 50Hz
- Timer : 1 to 15 min adjustable
- Dimension : 245 x 215 x 95mm
- Weight : 2.7 kg

 New No.8, Old No. J71, Karikalan Street
Jafferkhanpet, Chennai - 600 083

 +91-44-4858 2421

 +91-73586 20147 / 73586 20247

 alphameditronics@gmail.com

 www.alphameditronics.com
www.alphameditronics.in



ALPHA

ALPHA MEDITRONICS

PHYSIO & MEDICAL EQUIPMENTS MANUFACTURER & SALES



Alpha Meditronics



AlphaMeditronic



+91-73586 20147
+91-73586 20247



alphameditronics